


Calculating Molar Excess for Disposal



Example: 100 mL of a 0.1 M solution of potassium dichromate needs to be neutralized by a twofold molar excess of 4% sodium thiosulfate.

First, find the moles of potassium chromate solution:

$$=100 \text{ mL} \times 0.1 \text{ mol/L} \times 1 \text{ L} / 1000 \text{ mL} = 0.01 \text{ mol K}_2\text{Cr}_2\text{O}_7$$

Next, find the 4% sodium thiosulfate solution needs to be converted to concentration:

$$\% \text{ m/v} = \frac{\text{g of solute}}{\text{mL of solution}}$$

$$4\% \text{ m/v} = 4.0 \text{ g} / 100 \text{ mL}$$

$$\text{Concentration of neutralizing substance} = 4.0 \text{ g} / 100 \text{ mL} \times 1000 \text{ mL} / 1 \text{ L} \times \text{mol} / 158.11 \text{ g} = 0.25 \text{ mol/L}$$

Lastly, use the moles of solution and concentration of neutralization substance to find volume needed for reaction:

$$= 0.01 \text{ mol K}_2\text{Cr}_2\text{O}_7 \times 2 \text{ mol excess Na}_2\text{S}_2\text{O}_3 / 1 \text{ mol K}_2\text{Cr}_2\text{O}_7 \times 0.25 \text{ mol/L} \times 1000 \text{ mL} / 1 \text{ L} = 79 \text{ mL Na}_2\text{S}_2\text{O}_3$$

Conclusion: Add 79 mL of sodium thiosulfate to 100 mL of 0.1 M potassium dichromate to neutralize