

Data Table

Unknown	Trial	pH	pH (average)	p <i>K</i> _a	Unknown Identity
A	Sample #1				
	Sample #2				
B	Sample #1				
	Sample #2				
C	Sample #1				
	Sample #2				
D	Sample #1				
	Sample #2				
E	Sample #1				
	Sample #2				

Post-Lab Questions

- Average the pH readings for each trial (samples #1 and #2) to calculate the average p*K*_a value for the unknown weak acids and enter answers in the Data Table.
- Comment on the precision (reproducibility) of the p*K*_a determinations. Describe sources of experimental error and their likely effect on the measured p*K*_a (pH) values.
- The following table lists the identities of the unknowns in this experiment. Complete the table by calculating the p*K*_a value for each acid. *Note:* p*K*_a = -log*K*_a.

Acid	Formula	<i>K</i> _a	p <i>K</i> _a
Potassium dihydrogen phosphate	KH ₂ PO ₄	<i>K</i> _{a2} of H ₃ PO ₄ = 6.2 × 10 ⁻⁸	
Potassium hydrogen sulfate	KHSO ₄	<i>K</i> _{a2} of H ₂ SO ₄ = 1.0 × 10 ⁻²	
Potassium hydrogen phthalate	KHC ₈ H ₄ O ₄	<i>K</i> _{a2} of H ₂ C ₈ H ₄ O ₄ = 3.9 × 10 ⁻⁶	
Potassium hydrogen tartrate	KHC ₄ H ₄ O ₆	<i>K</i> _{a2} of H ₂ C ₄ H ₄ O ₆ = 4.6 × 10 ⁻⁵	
Acetylsalicylic acid	2-CH ₃ CO ₂ C ₆ H ₄ COOH	<i>K</i> _a = 3.2 × 10 ⁻⁴	

4. Compare the experimental pK_a value for each unknown with the literature values reported in Question 3. Determine the probable identity of each unknown and enter the answers in the Data Table.

5. Write separate equations for each unknown potassium salt dissolving in water and for the ionization reaction of the weak acid anion that each of these salts contains. (See Equations 7 and 8.)

6. Why was it not necessary to know the exact mass of each acid sample?

7. Why was it not necessary to know the exact concentration of the sodium hydroxide solution?

8. Why was it necessary to measure the exact volume of distilled water used to dissolve the acid, as well as the exact volume of solution transferred from the beaker to the Erlenmeyer flask?