

# Ruby Red Colloidal Gold

Copper, Silver and Gold Redox Reactions



## Introduction

From nanotech fibers and nanosensors to nanobots, nanotechnology has created so much “buzz” that it’s hard to tell where the science ends and the science fiction begins. Wherever it may lead in the future, the science of nanotechnology begins with solid particles called nanoparticles that are 1–100 nm in size. Shrinking the size of solid-phase particles to the nanometer scale—one billionth of a meter—changes their physical and chemical properties. The surprising properties of “colloidal gold” are a good example of this phenomenon. Whereas normal or “bulk” gold is bright, shiny, metallic yellow, colloidal gold nanoparticles are red or blue, and not at all shiny. Let’s investigate the preparation, properties, and uses of colloidal gold.

## Concepts

- Nanotechnology
- Colloids vs. solutions
- Redox reaction
- Metric measurements

## Materials

Hydrogen tetrachloroaurate (gold chloride) solution,  $\text{HAuCl}_4$ , 1 mM ( $1 \times 10^{-3}$  M), 20 mL

Sodium chloride solution,  $\text{NaCl}$ , 1 M, 5 mL

(Tri)sodium citrate solution,  $\text{Na}_3\text{C}_6\text{H}_5\text{O}_7$ , 1%, 2 mL

Water, distilled, 200 mL

Beaker, 100-mL

Beral-type pipets, graduated, 2

Ceramic fiber square (optional)

Cuvet or test tube, 13 × 100 mm

Erlenmeyer flask, 250- or 500-mL

Glass stirring rod

Graduated cylinder, 25-mL

Hot plate

Laser pointer or flashlight

Spectrophotometer or colorimeter (optional)

Test tubes, medium, 2

## Safety Precautions

*Dilute hydrogen tetrachloroaurate solution may be irritating to the eyes, skin, and the gastrointestinal tract. The potential health effects of nanoparticles have not been fully identified. Avoid contact of all chemicals with eyes and skin. Wear chemical splash goggles, chemical-resistant gloves, and a chemical-resistant apron. Please review current Material Safety Data Sheets for additional safety, handling, and disposal information.*

## Procedure

1. Measure 20 mL of 1 mM hydrogen tetrachloroaurate solution in a graduated cylinder, and pour the solution into a 250-mL Erlenmeyer flask.
2. Add distilled water to the 200-mL mark on the Erlenmeyer flask, and place a glass stirring rod into the solution.
3. Place the Erlenmeyer flask on a hot plate at a medium-high setting and heat the diluted gold chloride solution to gentle boiling.
4. Using a graduated, Beral-type pipet, add 2 mL of 1% trisodium citrate solution to the boiling solution in the Erlenmeyer flask.
5. Observe the color changes in the solution as gold(III) ions react with citrate. *(The pale yellow solution will gradually turn blue, violet, and then magenta or ruby-red.)*
6. Continue heating the solution at a gentle boil for approximately 10 minutes until the solution is ruby- or wine-red and the color no longer changes.

7. Remove the Erlenmeyer flask from the hot plate. Place the flask on a heat-resistant surface or ceramic fiber square and allow the solution to cool. Add distilled water to bring the total volume of liquid back up to 200 mL, if necessary.
8. When the solution has cooled to room temperature, pour some of the “colloidal gold” into a small beaker and observe its properties. (*The colloidal gold “solution” is an intense wine-red or reddish violet color. There are no solid-phase particles visible, but the liquid appears very slightly translucent rather than transparent.*)
9. Shine a laser pointer or flashlight through the colloidal gold solution and observe the “path” of the light through the solution. (*The beam of light becomes visible as it passes through the solution. This property, called the Tyndall effect, is due to light scattering.*)
10. Discuss what happens when light is shone through a true (transparent) solution such as water. (*A true solution is transparent—light passes through the solution and is projected on the other side, but the path of light in the solution itself is not visible.*)
11. Compare and contrast the following properties of a colloid versus a true solution: particle size, settling behavior, filtration, light scattering. (See the *Discussion* section.)
12. Pour some of the colloidal gold into two medium-size test tubes, filling each test tube about one-third full.
13. Add an equal volume of distilled water to the first test tube, and an equal volume of 1 M sodium chloride solution to the second test tube. Carefully swirl each test tube to mix the contents. Observe any color changes of the colloidal gold. (*The original red color of the colloidal gold changes to blue when it is diluted with sodium chloride solution. The red color of the colloidal gold simply fades to paler shade of red when diluted with water.*)
14. (*Optional*) Fill a cuvet approximately two-thirds full with the colloidal gold solution and measure the absorbance every 10 nm from 400 to 700 nm using a spectrophotometer. (See the *Discussion* section.)

## Disposal

Please consult your current *Flinn Scientific Catalog/Reference Manual* for general guidelines and specific procedures governing the disposal of laboratory waste. The colloidal gold solution is very stable and may be stored indefinitely. Keep the solution in a dark bottle to avoid exposure to light. Because of the unknown potential health hazards of colloidal gold, we do not recommend disposing of colloidal gold down the drain. The colloid may be broken by adding 6 M hydrochloric acid, which precipitates the gold. Solid gold may be disposed of in the trash according to Flinn Suggested Disposal Method #26a. Excess hydrogen tetrachloroaurate solution should be stored for future use.

## Tips

- Using a more concentrated solution of  $\text{HAuCl}_4$  in the colloidal gold preparation produces a cloudy, dark blue dispersion of gold nanoparticles. The blue mixture consists of larger particles (ca. 100 nm)—the visible absorbance shifts to longer wavelength, and there is considerable light scattering.
- The scanning tunneling microscope (STM) is the most important tool for studying the size and shape of nanoparticles. Using this non-optical, scanning probe microscope, scientists are able to “see” individual atoms and molecules at a resolution of 0.2 nm ( $2 \times 10^{-10}$  m). A tiny electrical probe or stylus is moved across a surface, producing a weak electrical current between the tip and the surface. The locations of atoms on the surface are visualized as regions of high electron density due to changes in the magnitude of the current or the position of the stylus. Heinrich Rohrer and Gerd Binnig of the IBM Research Laboratory in Zurich, Switzerland, received the Nobel Prize in Physics in 1986 for their invention of the STM. Visit the official Web site of the Nobel Foundation at [http://nobelprize.org/education\\_games/physics/microscopes/scanning/gallery/index.HTML](http://nobelprize.org/education_games/physics/microscopes/scanning/gallery/index.HTML) (accessed November 2006) to view an impressive gallery of “atomic” photos obtained with the STM.

## Discussion

Nanoscience or nanotechnology involves the preparation, characterization, and uses of nano-sized particles having dimensions in the 1–100 nm range. Nanoparticles have unique physical and chemical properties that are significantly different from the macroscopic properties of traditional or bulk solids. Many of these properties have taken on special importance in recent years as the applications of nanotechnology have intensively studied. In particular, the electronic, magnetic, and optical properties of nanoparticles have proven to be very useful in the creation of new products using nanotechnology. Quantum dots, for example, are nanocrystalline fluorescent semiconductors that are used in high definition DVD players and video game consoles.


Gold nanoparticles are one of the most widely used materials in nanotechnology, and they are certainly the oldest. Colloidal gold consists of gold(0) nanoparticles that range in size from approximately 5–50 nm and are uniformly dispersed in water. Michael Faraday (1791–1867) published a scientific report of the preparation and properties of colloidal gold in 1857, but the art of using colloidal gold as a decorative pigment in glassmaking dates back more than 1000 years. The Lycurgus Cup at the British Museum in London, England, a Roman artifact from the fourth century A.D., is the most famous example of ancient “gold ruby glass.” The presence of tiny, 70-nm crystals of silver and gold “dissolved” in the glass gives the cup a lustrous red appearance when light shines through it. This demonstrates the most striking and beautiful feature of gold nanoparticles—their color. Depending on the size and shape of the particles, the color of gold nanoparticles varies from red to purple. The optical properties of gold nanoparticles are not only unique, they are also useful, providing the basis for commercial products such as medical diagnostic kits for HIV detection, biosensors for DNA analysis, lasers, and optical filters. From a chemical standpoint, gold nanoparticles have a very active surface chemistry and are thus valuable catalysts for pollution control, fuel cells, and the synthesis of specialty chemicals. The ability of various electron-donating groups to bind to the surface of gold nanoparticles is being investigated as a way to direct the self-assembly of complex structures for futuristic nano-computers and other nano-electronic devices.

The most common method for the preparation of colloidal gold involves the reduction of gold(III) ions by citrate ions in dilute (0.1 mM) aqueous solution. The gold(III) ions are usually added to water in the form of hydrogen tetrachloroaurate (HAuCl<sub>4</sub>). Citrate ions (C<sub>6</sub>H<sub>5</sub>O<sub>7</sub><sup>3-</sup>), which are most likely oxidized to acetone dicarboxylate ions (C<sub>5</sub>H<sub>4</sub>O<sub>4</sub><sup>2-</sup>) in the process, act as a two-electron reducing agent. The half reactions for the synthesis of colloidal gold are given in Equations 1 and 2.



The average diameter of gold nanoparticles produced by this method depends on temperature and on the concentration ratio of gold(III) ions and citrate ions in solution. The gold nanoparticles are stabilized by the presence of citrate ions adsorbed on the surface of the particles. Adsorption of citrate ions gives the gold particles an overall negative charge and is the principal factor responsible for the formation of a stable colloid. Mutual repulsion of the small, negatively charged particles prevents them from coagulating to form larger particles that might eventually settle out of solution.

The visible absorption spectrum of the colloidal gold produced in this demonstration is shown in Figure 1. The wavelength of maximum absorbance at 520 nm correlates with the formation of gold nanoparticles having an average diameter of 20–40 nm. The “peak width at half-maximum” for this colloidal gold preparation is quite broad (130 nm) and indicative of a fairly wide distribution of particle sizes around the mean. In general, the wavelength of maximum absorbance shifts to higher wavelength (>520 nm) when the mean particle size increases above 40 nm, and the peak width increases when there is a larger variation in particle sizes. The color change of the gold colloid from red to blue when sodium chloride is added illustrates this effect. Adding NaCl, a strong electrolyte, shields the negative charges of the colloidal gold nanoparticles and causes them to clump together to form larger particles.



**Figure 1.** Visible Spectrum of Colloidal Gold.

The absorbance of visible light by gold and other metal nanoparticles has been attributed to a unique phenomenon called *surface plasmon resonance* (SPR). This phenomenon is very different from the “normal” visible spectrum of colored organic dye molecules, for example, which is due to the promotion of electrons from the ground state to an excited state when light of a specific wavelength is absorbed. SPR is defined as the “collective oscillation of conduction band electrons resulting from the interaction with electromagnetic radiation.” In laymen’s terms, the incoming electromagnetic radiation induces the formation of a dipole on the surface of the nanoparticle, which then oscillates in phase or in resonance with the electric field of the incoming light. This occurs at a specific frequency (and wavelength or color) of light, depending on the size, shape, and form of the nanoparticles. As noted above, dispersions containing small gold nanoparticles are red. When the nanoparticles begin to aggregate further, their color changes from red to blue and then to purple.

Solutions and colloids, which differ in the size of the particles that are dispersed throughout a continuous phase, are distinguished from one another primarily in terms of their properties. Colloids, for example, may be defined as mixtures in which the dispersed particles are small enough to pass through a filter but too large to pass through a semipermeable membrane. The particles in a colloid are large enough that they will reflect or scatter light in all directions. The scattering of light by particles in a mixture is called the Tyndall effect and makes it possible to view a beam of light as it passes through a colloid or a suspension. In a true solution the dispersed particles are too small to scatter visible light. The following table summarizes the properties of solutions, colloids, and suspensions. Notice that the particle size range for each type of mixture is just that, a range, and not an absolute or fixed value. There is thus a continuum of properties for solutions, colloids, and suspensions.

Property	Solution	Colloid	Suspension
Particle Size	0.1–1 nm (atoms, ions, small molecules)	1–200 nm	>200 nm (aggregates of large molecules)
Light Scattering	None	Tyndall effect	Tyndall effect
Settling Behavior	Stable, does not separate.	Stable, does not separate	Particles separate on standing.
Filtration	Particles pass through filter.	Particles pass through filter.	Particles do not pass through filter.

## Connecting to the National Standards

This laboratory activity relates to the following National Science Education Standards (1996):

**Unifying Concepts and Processes: Grades K–12**

Systems, order, and organization  
Evidence, models, and explanation

**Content Standards: Grades 5–8**

Content Standard B: Physical Science; properties and changes of properties in matter  
Content Standard F: Science in Personal and Social Perspectives; science and technology in society

**Content Standards: Grades 9–12**

Content Standard B: Physical Science; structure of atoms; structure and properties of matter; interactions of energy and matter  
Content Standard F: Science in Personal and Social Perspectives; science and technology in local, national, and global challenges

## Answers to Worksheet Questions *(Student answers will vary.)*

- The line drawn on the side of this page is 220 mm long. Assume that this line represents the average diameter of a red blood cell, which is 1 micrometer or micron ( $1 \mu\text{m} = 1 \times 10^{-6} \text{ m}$ ). (a) What is the diameter of a red blood cell in *nanometers*? (b) Calculate the length of the line segment in *millimeters* that would be needed to represent an object that is one *nanometer* in diameter on this line.

a.  $\frac{1 \times 10^9 \text{ nm}}{1 \text{ m}} = 1 \times 10^{-6} \text{ m} \times \frac{220 \text{ mm}}{1 \text{ m}} = 220 \times 10^3 \text{ nm}$  (The diameter of a red blood cell is 1000 nm.)

b. If 1000 nm corresponds to 220 mm on the line, then  $220/1000 = 0.22 \text{ mm}$  would be needed to represent an object that is 1 nm in diameter.

- Mark off line segments of the appropriate length on the line to represent (a) the average size of the influenza virus (100 nm), (b) a colloidal gold nanoparticle (40 nm), and (c) the width of the DNA double helix (2 nm).

*The length of the line segments should be (a) 22 mm for the influenza virus, (b) 9 mm for the gold nanoparticle, and (c) 0.4 mm for the DNA helix.*

- The colloidal gold in this demonstration was prepared starting with 20 mL of  $1 \times 10^{-3} \text{ M H AuCl}_4$ . (a) How many grams of gold (Au = 197 g/mole) are contained in the flask of colloidal gold? (b) At a current price of \$580 per troy ounce (1 troy ounce = 31.1 g) for gold, how much is the gold in the flask worth?

a. *The mole ratio is one mole of gold produced per mole of  $\text{H AuCl}_4$  added to the flask.*

$$1 \times 10^{-3} \text{ moles/L} \times 0.020 \text{ L} \times 197 \text{ g/mole} = 0.0039 \text{ g of gold}$$

b.  $0.0039 \text{ g} \times (1 \text{ Troy ounce}/31.1 \text{ g}) \times (\$580/\text{Troy ounce}) = \$0.073$  (The gold is worth seven cents!)

- Estimate the number of gold atoms in a single gold nanoparticle that is 40 nm in diameter. Use the following assumptions: (a) The atomic radius for gold is 0.15 nm. (b) Both particles are the shape of a sphere. The volume of a sphere is  $\frac{4}{3}\pi r^3$ , where r is the radius. (c) Only 74% of the total volume of the nanoparticle is physically occupied by gold atoms. (The rest of the volume is “empty space” between atoms. 74% is the maximum “packing efficiency” of spheres in any crystal lattice.)

*The total volume of the gold nanoparticle is equal to  $\frac{4}{3}(3.14)(20 \text{ nm})^3$  or  $3.3 \times 10^4 \text{ nm}^3$*

*The volume of a gold atom is equal to  $\frac{4}{3}(3.14)(0.15 \text{ nm})^3$  or  $0.014 \text{ nm}^3$*

*The effective volume of the nanoparticle that is occupied by gold atoms is 74% of the total volume, or  $(0.74)(3.3 \times 10^4 \text{ nm}^3) = 2.4 \times 10^4 \text{ nm}^3$ .*

*The approximate number of gold atoms can be obtained by dividing the effective volume of the nanoparticle by the volume of one gold atom.*

$$\frac{2.4 \times 10^4 \text{ nm}^3}{0.014 \text{ nm}^3/\text{gold atom}} = 1.7 \times 10^6 \text{ gold atoms}$$

*There are almost two million gold atoms in a single nanoparticle!*

### Reference

Liz-Marzán, Luis M. “Nanometals: Formation and Color” *Materials Today* **2004**, 7, 26–31. May be downloaded from the author’s university web site at [webs.uvigo.es/coloides/nano/pdf/MT-04.pdf](http://webs.uvigo.es/coloides/nano/pdf/MT-04.pdf) (accessed December 2006).

Turkevich, J. “Colloidal Gold. Part I. Historical and Preparative Aspects, Morphology and Structure” *Gold Bulletin* **1985**, 18, 86–91. May be downloaded from the Gold Bulletin Web site at [http://www.goldbulletin.org/downloads/Turkevich\\_3\\_18.pdf](http://www.goldbulletin.org/downloads/Turkevich_3_18.pdf) (accessed December 2006).

### Flinn Scientific—Teaching Chemistry™ eLearning Video Series

A video of the *Ruby Red Colloidal Gold* activity, presented by Irene Cesa, is available in *Copper, Silver and Gold Redox Reactions* and in *Industrial Chemistry*, part of the Flinn Scientific—Teaching Chemistry eLearning Video Series.

### Materials for *Ruby Red Colloidal Gold* are available from Flinn Scientific, Inc.

Materials required to perform this activity are available in the *Ruby-Red Colloidal Gold—Chemical Demonstration Kit* available from Flinn Scientific.

Catalog No.	Description
AP7117	Ruby-Red Colloidal Gold—Chemical Demonstration Kit

Consult your *Flinn Scientific Catalog/Reference Manual* for current prices.

# Ruby Red Colloidal Gold Worksheet

1. The line drawn on the side of this page is 220 mm long. Assume that this line represents the average diameter of a red blood cell, which is 1 micrometer or micron ( $1 \mu\text{m} = 1 \times 10^{-6} \text{ m}$ ). (a) What is the diameter of a red blood cell in *nanometers* ( $1 \text{ nm} = 1 \times 10^{-9} \text{ m}$ )? (b) Calculate the length of the line segment in *millimeters* that would be needed to represent an object that is one *nanometer* in diameter on this line.
  
2. Mark off line segments of the appropriate length on the line to represent (a) the average size of the influenza virus (100 nm), (b) a colloidal gold nanoparticle (40 nm), and (c) the width of the DNA double helix (2 nm).
  
3. The colloidal gold in this demonstration was prepared starting with 20 mL of  $1 \times 10^{-3} \text{ M H AuCl}_4$ . (a) How many grams of gold ( $\text{Au} = 197 \text{ g/mole}$ ) are contained in the flask of colloidal gold? (b) At a current price of \$580 per troy ounce (1 troy ounce = 31.1 g) for gold, how much is the gold in the flask worth?
  
4. *Estimate* the number of gold atoms in a single gold nanoparticle that is 40 nm in diameter. Use the following assumptions: (a) The atomic radius for gold is 0.15 nm. (b) Both particles are the shape of a sphere. The volume of a sphere is  $\frac{4}{3}\pi r^3$ , where  $r$  is the radius. (c) Only 74% of the total volume of the nanoparticle is physically occupied by gold atoms. (The rest of the volume is “empty space” between atoms. 74% is the maximum “packing efficiency” of spheres in any crystal lattice.)