

Measuring with Laser Light Worksheet

Data Table(s) or Observations

Object	Distance from object to screen, L (mm)	Distance between dark spaces, Δy (mm)	Wavelength of laser light, λ (mm)	Object width (mm)	Object width (μm)
Fishing line					
Copper wire					
Human hair					

Post-Lab Calculations *(Show your work on a separate sheet of paper.)*

- For each object, measure the distance in mm between pieces of tape (see Figure 7) on the paper screen. Record the results in the data table.
- Record the wavelength of the laser light in mm. *Note:* Most lasers have a label with the wavelength recorded in nm. Convert nm to mm. If your laser pointer does not indicate the wavelength, check with your instructor.
- Calculate the width of each object in mm using Equation 1. Round to the nearest hundredth of a mm and record in the data table.
- Convert the width of each object from mm to μm and record the results.

Post-Lab Questions and Analysis

- List the measured objects from the data table in order from smallest width to largest. How did the diffraction pattern change from one object to the next?
- The diameter of human hair varies, but is usually in the range of 20–180 μm . Did the experimental value obtained fall within this range? Compare your results with other groups. Does there seem to be a relationship between hair color and width? Explain.
- The diameter of a 30-gauge copper wire is 0.255 mm. a) How does the measured width of the copper wire compare to the accepted width? b) Use Equation 2 to calculate the percent error between the measured and accepted values for the width of the wire. c) What are some possible sources of error in this experiment?

$$\text{Percent error} = \frac{|\text{measured value} - \text{accepted value}|}{\text{accepted value}} \times 100\% \quad \text{Equation 2}$$