

Empirical Formula of Copper Carbonate Worksheet

Data Table. Part I. Gas Evolution Method

| | |
|---|--|
| Mass of clean and dry Erlenmeyer flask | |
| Mass of flask and basic copper carbonate sample | |
| Mass of basic copper carbonate analyzed | |
| Mass of cylinder and HCl | |
| Mass of cylinder after HCl was added to reaction flask | |
| Mass of HCl added to flask | |
| Mass of flask + basic copper carbonate sample + HCl | |
| Mass of flask + final solution after CO ₂ loss | |
| Mass of released CO ₂ | |
| Percent CO ₂ in basic copper carbonate sample | |

Data Table. Part II. Colorimetric Comparison

| Test Tube | 1 | 2 | 3 | 4 | 5 | 6 | Unknown |
|--|---------|--------|--------|--------|--------|--------|---|
| Volume of CuSO ₄ stock solution | 10.0 mL | 8.0 mL | 7.0 mL | 6.0 mL | 4.0 mL | 2.0 mL | 1.30 g of basic copper carbonate diluted to 100.0 mL with 0.2 M sulfuric acid |
| Final volume of diluted solution (mL) | 10.0 | 10.0 | 10.0 | 10.0 | 10.0 | 10.0 | 10.0 |
| Color comparison (rank solutions from lightest blue = 1 to deepest blue = 6) | 6 | 5 | 4 | 3 | 2 | 1 | The known solution whose color matches the unknown = |
| Copper ion concentration in standard solution | 0.20 M | 0.16 M | 0.14 M | 0.12 M | 0.08 M | 0.04 M | Estimated concentration of copper ion in the unknown solution = |

Post-Lab Analysis and Calculations

1. Show the calculations for the gas evolution method Data Table.

- a. Mass of basic copper carbonate analyzed
- b. Mass of HCl added to flask
- c. Mass of flask + basic copper carbonate sample + HCl
- d. Mass of released CO₂
- e. Percent CO₂ in basic copper carbonate sample

2. Calculate the percent copper based on the results from Part II, the colorimetric comparison. Example: if the color of the unknown basic copper carbonate solution most closely resembled test tube #2, which has a concentration of 0.16 M, the calculation would be:

$$\frac{0.16 \text{ moles}}{\text{L}} \times 0.1 \text{ L} \times \frac{63.55 \text{ g Cu}}{\text{mol}} = 1.02 \text{ g Cu in the 100 mL unknown solution}$$

Unknown basic copper carbonate solution size is 100 mL.

$$\frac{1.02 \text{ g Cu}}{1.30 \text{ g copper carbonate}} \times 100\% = 78\% \text{ Cu}$$

3. Fill in the table below and identify the most likely form of basic copper carbonate that was tested in this lab. Write a short paragraph explaining your choice and describe the supporting evidence quantitatively.

| | % CO ₂ theoretical | % CO ₂ experimental | % copper theoretical | % copper experimental |
|---|-------------------------------|--------------------------------|----------------------|-----------------------|
| Cu ₂ (OH) ₂ CO ₃ | | | | |
| Cu ₃ (OH) ₂ (CO ₃) ₂ | | | | |