



Online Chemventory™ Tutorials

Adding SDS links for non-Flinn chemicals


Any chemical added to your inventory list from the following three sources will include a link to a Flinn safety data sheet (SDS):

- Chemicals added to your inventory via results from a keyword search.
- Chemicals added to your inventory via the Flinn chemical checkbox selection list.
- Chemicals (*with Flinn catalog numbers listed*) added to your inventory via an import from Chemventory 8.0 or 5.5 disc based software.

Chemicals added manually will not have an SDS automatically associated with them. You can add a link to an SDS for your chemical by following the steps below.

Note: Flinn Online Chemventory™ will allow you to include a link to any safety data sheet in HTML or PDF format that resides on an outside website. Flinn Scientific is not responsible for the content, updating the content, or maintaining broken links for any non-Flinn SDS linked to your Online Chemventory™ database.

- Type the chemical name and the words “safety data sheet” into your search engine (Google®, Yahoo®, Bing®). In the example below we have entered “Ammonium Nitrate Safety Data Sheet” in a Google® search. If you have the company name for the chemical SDS you are looking to add, for example Sigma Aldrich, you may also enter that name to provide better search results.



Continued on next page



- 1) Your top 10 search results will most likely provide at least 5 safety data sheet PDF options. Most search engines will include a notation next to the page title if that link is a PDF as shown in the screen shot below. Click on the link for the PDF to open that page.

[PDF] [Ammonium nitrate MSDS - ScienceLab](#)
www.sciencelab.com/msds.php?msdsid=9927336/ ▾
Oct 10, 2005 - Material **Safety Data Sheet. Ammonium nitrate** MSDS. Section 1:
Chemical Product and Company Identification. Product Name: Ammonium ...

[Ammonium Nitrate](#)
avogadro.chem.iastate.edu/msds/nh4no3.htm ▾ Iowa State University ▾
Material **Safety Data Sheet ...** MSDS Name: **Ammonium Nitrate** Catalog ...
Synonyms: Nitric acid ammonium salt, nitram, ammonium salt, norway saltpeter

[PDF] [Safety Data Sheet - Orica MSDS](#)
msds.orica.com/pdf/shess-en-cds-010-000032017701.pdf ▾ Orica ▾
Safety Data Sheet. 1. IDENTIFICATION OF THE MATERIAL AND SUPPLIER.
Product Name: **AMMONIUM NITRATE**. Other name(s): Nitric acid ammonium salt; ...

[PDF] [Material Safety Data Sheet Calcium Ammonium Nitrate](#)
<http://www.lakeland.edu/.../MSDS/.../Instant%20MSDS.pdf> ▾ Lakeland College ▾
Material **Safety Data Sheet**. Product Code Number: 31246, 31246-b, 31359, 31346, 91640. Issue Date: December 2010.

[PDF] [Ammonium Nitrate - CSBP](#)
www.csbp.com.au/.../ammonium-nitrate/.../1210898960_ammonium_nitr... ▾
May 16, 2008 - MATERIAL **SAFETY DATA SHEET. Ammonium Nitrate**. Section 1 –
Identification of the Material and Supplier. Product Name: **Ammonium nitrate**.

[PDF] [SC "Achema" Safety data sheet Ammonium nitrate](#)
www.achema.com/.../Ammonium%20nitrate%20N34.../2011.12.30%20... ▾
Dec 30, 2011 - Use of **ammonium nitrate** as an intermediate to synthesize other ...
Person responsible for the **Safety Data Sheet** (with e-mail address): Vidas ...

[PDF] [Ammonium Nitrate Amtrate MSDS - CF Industries, Inc.](#)
www.cfindustries.com/.../Ammonium-Nitrate... ▾ CF Industries Holdings, Inc. ▾
Page 1 of 9. Material **Safety Data Sheet. Ammonium Nitrate**. MSDS Number 004
(Revised April 23, 2013). 9 Pages. 1. CHEMICAL PRODUCT and EMERGENCY ...

[PDF] [Material Safety Data Sheet: AMMONIUM NITRATE ... - CALAMCO](#)
www.calamco.com/pdf/ammonium_nitrate_msd.pdf ▾
Material **Safety Data Sheet: AMMONIUM NITRATE SOLUTION (20-0-0 AN)**. Last
Revision 01/02/13 Page 1 of 3. SECTION 1: CHEMICAL PRODUCT ...

Continued on next page



- Highlight the entire URL once you have found a safety data sheet link that you would like to add to your Online Chemventory™ database. After the URL is highlighted, right click on it and select copy.

Note: The entire URL including the *http://* portion must be highlighted and copied.

Material Safety Data Sheet

Product Code Numbers: 715, 718-b, 31259, 31246, 31246-b, 31337, 31346, 91640

Issue Date: November 2010

Calcium Ammonium Nitrate
Conforms to ISO 11014-1 and 91/155/EC

1. Chemical Product and Company Identification

Chemical product name: Calcium ammonium nitrate
Common chemical name:
Synonyms: CAN/Nutramon
Supplier's name:
Supplier's address:

2. Composition/Information on Ingredients

Name	CAS number	% by weight	Symbol	R phrases
Ammonium Nitrate	6484-52-2	max. 78	—	—

Additional Information
Contains:

Name	CAS number
Magnesium Nitrate	13446-18-9

Continued on next page



- Return to the Online Chemventory™ program and follow the next four steps to insert the URL hyperlink that you copied in the previous step:

1) Click the name of the chemical for which you would like to add the safety data sheet URL link.

2) From the Chemical Information page, click the “Edit Info” button in the gray navigation bar.

Second Floor Store Room chemventory@test.com Sign Up Log Out

FLINN SCIENTIFIC, INC. "Your Safer Source for Science Supplies" **Online Chemventory™**

Home Inventory Purchase Lists Label Print Queue

+ Add Chemical + Manage Users + Import + Export Need Help?

Inventory

Name: Second Floor Store Room 4 Chemicals Currently Listed
School: Flinn Scientific Last Updated: Feb 17, 2014
ID: 10011

Email Selection Print Selection Remove Selection

All	Chemical Name	School	Store Room	Shelf	Amount	SDS
<input type="checkbox"/>	Aluminum, Granular (5 Mesh)	Flinn Sci.	Second Floor	I5 top shelf	1000 mL	
<input type="checkbox"/>	Hydrochloric Acid (12 M)	Flinn	Upstairs Chem.	Acid Cabinet	800 mL	
<input type="checkbox"/>	Adenosine Triphosphate, Disodium	Flinn	Upstairs chem	Right H2	2 g	
<input type="checkbox"/>	Dextrose	Flinn	upstairs Chem	Right H4	100 g	

Second Floor Store Room chemventory@test.com Sign Up Log Out

FLINN SCIENTIFIC, INC. "Your Safer Source for Science Supplies" **Online Chemventory™**

Home Inventory Purchase Lists Label Print Queue

Edit Info Add to Label Print Queue Add to Purchase List Print Chemical Info

Dextrose Agar

Second Floor Store Room

Chemical Information

As of Feb 17, 2014

Name:	Dextrose Agar	Unit of Measure:	g
Formula:		Minimum Reorder Amount:	g
Synonym:		GHS Signal Word:	
Company:	Flinn Scientific	Hazard Statement:	
Catalog #:	D0049	This chemical is considered nonhazardous according to GHS classifications for the Hazard Communication Standard. Treat all laboratory chemicals with caution.	
Grade:	NONE	GHS Pictograms:	
Disposal #:	26a	SDS Information:	http://www.flinnsci.com/Documents/SDS/D/DextroseAgar.pdf
CAS #:	None Established	Notes:	
Class:	ORG		
Compatible Family:	O M		

Inventory Information

School:	Flinn	Amount:	100 g
Store Room:	upstairs Chem		
Shelf:	Right H4		


Continued on next page



3) Locate the SDS Information Link data field on the Chemical Information Edit Info page.

4) Right click inside this data field and select "Paste." After the URL link has been added, scroll to the bottom of the page and click Submit.

The screenshot shows the 'Cold Pack, Instant' Chemical Information page. The 'SDS Information Link' field is circled in red. A red arrow points from this field to a zoomed-in view of the field and its context menu.



Continued on next page



You will now be able to click the SDS Icon link in your inventory list to open the safety data sheet that was linked.

Second Floor Store Room chemventory@test.com Sign Up Log Out

FLINN SCIENTIFIC, INC. "Your Safer Source for Science Supplies" **Online Chemventory™**





Home Inventory Purchase Lists Label Print Queue

+ Add Chemical Manage Users Import Export Need Help?

Inventory

Name: Second Floor Store Room 4 Chemicals Currently Listed
School: Flinn Scientific Last Updated: Feb 17, 2014
ID: 10011

Email Selection Print Selection Remove Selection

All	Chemical Name	School	Store Room	Shelf	Amount	SDS
<input type="checkbox"/>	Aluminum, Granular (5 Mesh)	Flinn Sci.	Second Floor	15 top shelf	1000 mL	
<input type="checkbox"/>	Hydrochloric Acid (12 M)	Flinn	Upstairs Chem.	Acid Cabinet	800 mL	
<input type="checkbox"/>	Adenosine Triphosphate, Disodium Salt	Flinn	Upstairs chem	Right H2	2 g	
<input type="checkbox"/>	Dextrose Agar	Flinn	upstairs Chem	Right H4	100 g	

Material Safety Data Sheet

Product Code Numbers: 715, 716, 718, 718-b, 31259, 31246, 31246-b, 31359, 31346, 91640

Issue Date: November 2010

Calcium Ammonium Nitrate
Conforms to ISO 11014-1 and 91/155/EC

1. Chemical Product and Company Identification

Chemical product name: Calcium ammonium nitrate
Common chemical name:
Synonyms: CAN/Nutramon
Supplier's name:
Supplier's address:

2. Composition/Information on Ingredients

Name	CAS number	% by weight	Symbol	R phrases
Ammonium Nitrate	6484-52-2	max. 78	-	-

Additional Information

Contains:

Name	CAS number
Magnesium Nitrate	13446-18-9

Note: Flinn Scientific is not responsible for the content, updating the content, or maintaining broken links for any non-Flinn SDS linked to your Online Chemventory™ database. We recommend that you check with the original supplier of the chemical on an annual basis to see if updates have been made to any linked safety data sheet.