

10 Care & Maintenance Tips for Water Baths



1 Unplug when not in use and before cleaning. Allow it to cool before proceeding.

2 Clean with a mild detergent and soft cloth. Do not use bleach, as it will corrode the stainless steel.

3 Do not fill it with tap water or ultra-pure deionized water. Use distilled or deionized water.

4 Do not move the bath when it is filled with hot water. Do not touch the bath's hot surfaces when in use.

5 Do not operate with a damaged cord or plug. Do not let the cord touch hot surfaces.

6 Completely drain large baths using the drain outlet. Invert small baths carefully over a sink to drain.

7 Drain, clean and refill your water bath on a weekly basis if used regularly.

8 Keep water level at least half full. Ideally, water should cover the heater to avoid running it dry.

9 Do not place on or near a hot gas or electric burner, or on a heated surface.

10 Always turn off water when not in use. Don't leave on overnight.