



DANGER! Highly flammable liquid and vapor. Keep away from heat, sparks, and open flame. Toxic if swallowed, in contact with skin or inhaled. Causes damage to organs. Flash point 9.7 °C. TLV: 200 ppm.

FIRST AID: Flush affected area with water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF EXPOSED: Call a POISON CENTER or physician. IF SWALLOWED: Immediately call a POISON CENTER or physician.

FLINN
SCIENTIFIC

“Your Safer Source
for Science”

M0054

METHYL ALCOHOL

500 mL

absolute, reagent, CH_3OH , F.W. 32.04

★ **HAZARD ALERT:** Highly flammable liquid and vapor. Toxic if swallowed, in contact with skin, or if inhaled. Causes damage to organs.

770 N. Raddant Road, Batavia, IL 60510 U.S.A.
630-879-6900
© 2016 All Rights Reserved.

LOT:

STORAGE: Organic #2 inside a dedicated flammables cabinet. If a flammables cabinet is not available, store in a Flinn Safe-Stor Can.

ORGANIC #2



DISPOSAL: #18a

SHELF LIFE: Good if stored safely.

SOLUBLE: Miscible with water and alcohol.

CAS NO: 67-56-1

UN1230



NFPA

