# FLINN SCIENTIFIC, INC. Safety Data Sheet (SDS)

**SDS #:** 498.10

Revision Date: March 25, 2014

### SECTION 1 — CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

## Mercury(II) lodide Solution

Flinn Scientific, Inc. P.O. Box 219, Batavia, IL 60510 (800) 452-1261

CHEMTREC Emergency Phone Number: (800) 424-9300

Signal Word

**DANGER** 

## Pictograms

#### **SECTION 2 — HAZARDS IDENTIFICATION**

Hazard class: Acute toxicity, oral and inhalation (Category 2). Fatal if swallowed or inhaled (H300+H330). Do not eat, drink or smoke when using this product (P270). Do not breath mist, vapors or spray (P260).

Hazard class: Acute toxicity, dermal (Category 1). Fatal in contact with skin (H310).

Hazard class: Specific target organ toxicity - repeated exposure (Category 2). May cause damage to the kidneys, gastrointestinal tract, and nerves through prolonged or repeated exposure if inhaled, ingested, or absorbed through the skin (H373).

Hazard class: Hazardous to the aquatic environment, acute and chronic toxicity (Category 1). Very toxic to aquatic life with long lasting effects (H410). Avoid release to the environment (P273).





SECTION 3 — COMPOSITION, INFORMATION ON INGREDIENTS

Component Name	CAS Number	Formula	Formula Weight	Concentration
Mercury(II) iodide	7774-29-0	$\mathrm{HgI}_2$	454.43	4.5%
Water	7732-18-5	$H_2O$	18.00	95.5%

## SECTION 4 — FIRST AID MEASURES

Immediately call a POISON CENTER or physician (P310).

If inhaled: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do so. Continue rinsing.

**If on skin:** Gently wash with plenty of soap and water (P302 + P350).

If swallowed: Rinse mouth. Call a POISON CENTER or physician if you feel unwell.

#### **SECTION 5 — FIRE FIGHTING MEASURES**

Nonflammable, noncombustible solution.

**In case of fire:** Use a tri-class dry chemical fire extinguisher.

NFPA CODE

None established

#### **SECTION 6 — ACCIDENTAL RELEASE MEASURES**

Remove all ignition sources and water. Sweep up the spill, place in a sealed bag or container, and dispose. Ventilate area and wash spill site after material pickup is complete. See Sections 8 and 13 for further information.

**Safety Data Sheet** 

Mercury(II) Iodide Solution

**SDS #:** 498.10

Revision Date: March 25, 2014

#### SECTION 7 — HANDLING AND STORAGE

Flinn Suggested Chemical Storage Pattern: Inorganic #2. Store with acetates, halides, sulfates, sulfates, thiosulfates and phosphates. Store in a cool, dry place. Protect from light. Store in a Flinn Chem-Saf<sup>TM</sup> bag. Store in a poisons cabinet.

Keep container tightly closed (P233). Use only in a hood or well-ventilated area (P271).

#### SECTION 8 — EXPOSURE CONTROLS, PERSONAL PROTECTION

Wear protective gloves, protective clothing, and eye protection (P280). Wash hands thoroughly after handling (P264). Use only in a hood or well-ventilated area (P271).

Exposure guidelines: (as mercury iodide) PEL 0.1 mg/m<sup>3</sup> (OSHA); TLV 0.025 mg/m<sup>3</sup> (ACGIH) Readily absorbed through the skin.

#### - PHYSICAL AND CHEMICAL PROPERTIES **SECTION 9 -**

Red, tetragonal crystals. Odorless. Soluble: Boiling alcohol and solutions of sodium thiosulfate or

potassium iodide. Insoluble in water.

Boiling point: 354 °C Melting point: 259 °C Specific gravity: 6.28

#### **SECTION 10 — STABILITY AND REACTIVITY**

Avoid strong oxidizing and reducing agents.

Shelf life: Good, if stored properly.

#### SECTION 11 — TOXICOLOGICAL INFORMATION

Acute effects: Highly toxic, severe irritant ORL-RAT LD<sub>50</sub>: 18 mg/kg for mercury(II) iodide

Chronic effects: Neuro toxic effects, neurological hazard IHL-RAT LC<sub>50</sub>: N.A.

Target organs: Kidneys, GI system, nerves SKN-RAT LD<sub>50</sub>: 75 mg/kg for mercury(II) iodide

N.A. = Not available, not all health aspects of this substance have been fully investigated.

#### **SECTION 12 — ECOLOGICAL INFORMATION**

Very toxic to aquatic life with long lasting effects. Avoid release into the environment.

## <u>SECTION 13 — DISPOSAL CONSIDERATIONS</u>

Please review all federal, state and local regulations that may apply before proceeding.

Flinn Suggested Disposal Method #27f is one option.

#### **SECTION 14 — TRANSPORT INFORMATION**

Shipping name: Mercury Iodide; Hazard class: 6.1, Poison; UN number: UN1638

N/A = Not applicable

#### **SECTION 15 — REGULATORY INFORMATION**

TSCA-listed, EINECS-listed (231-873-8), RCRA code D009.

#### **SECTION 16 — OTHER INFORMATION**

This Safety Data Sheet (SDS) is for guidance and is based upon information and tests believed to be reliable. Flinn Scientific, Inc. makes no guarantee of the accuracy or completeness of the data and shall not be liable for any damages relating thereto. The data is offered solely for your consideration, investigation, and verification. The data should not be confused with local, state, federal or insurance mandates, regulations, or requirements and CONSTITUTE NO WARRANTY. Any use of this data and information must be determined by the science instructor to be in accordance with applicable local, state or federal laws and regulations. The conditions or methods of handling, storage, use and disposal of the product(s) described are beyond the control of Flinn Scientific, Inc. and may be beyond our knowledge. FOR THIS AND OTHER REASONS, WE DO NOT ASSUME RESPONSIBILITY AND EXPRESSLY DISCLAIM LIABILITY FOR LOSS, DAMAGE OR EXPENSE ARISING OUT OF OR IN ANY WAY CONNECTED WITH THE HANDLING, STORAGE, USE OR DISPOSAL OF THIS PRODUCT(S)

Consult your copy of the Flinn Science Catalog/Reference Manual for additional information about laboratory chemicals.

Revision Date: March 25, 2014