# FLINN SCIENTIFIC, INC. Safety Data Sheet (SDS)

**SDS #:** 505.10

**DANGER** 

Revision Date: March 25, 2014

#### SECTION 1 — CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

# **Mercury(I) Nitrate Solution**

Flinn Scientific, Inc. P.O. Box 219, Batavia, IL 60510 (800) 452-1261

CHEMTREC Emergency Phone Number: (800) 424-9300 Signal Word

**Pictograms** 







#### **SECTION 2 — HAZARDS IDENTIFICATION**

Hazard class: Acute toxicity, oral, inhalation, and dermal (Category 2). Fatal if swallowed, inhaled, or in contact with skin (H300+H310+H330). Do not eat, drink or smoke when using this product (P270). Do not breath mist, vapors or spray (P260).

Hazard class: Specific target organ toxicity - repeated exposure (Category 2). May cause damage to the kidneys, nerves, gastrointestinal tract, and blood through prolonged or repeated exposure if inhaled, ingested, or absorbed through the skin (H373).

Hazard class: Oxidizing liquids; Oxidizing solids (Category 2, 3). May intensify fire; oxidizer (H272). Keep away from heat, sparks, open flames, and hot surfaces. No smoking (P210).

Hazard class: Skin corrosion or irritation (Category 1). Causes severe skin burns and eye damage (H314). Hazard class: Hazardous to the aquatic environment, acute and chronic toxicity (Category 1). Very toxic to aquatic life with long lasting effects (H410). Avoid release to the environment (P273).

Avoid contact with acetic acid and readily oxidized substances.

Industrial exposure to nitric acid vapors and mists is listed as a known human carcinogen by IARC (IARC-1).

SECTION 3 — COMPOSITION, INFORMATION ON INGREDIENTS

Component Name	CAS Number	Formula	Formula Weight	Concentration
Mercury(I) nitrate	7782-86-7	HgNO <sub>3</sub>	525.19	24%
Nitric acid	7697-37-2	HNO <sub>3</sub>	63.00	24%
Water	7732-18-5	H <sub>2</sub> O	18.00	52%

## SECTION 4 — FIRST AID MEASURES

Immediately call a POISON CENTER or physician (P310).

If inhaled: Remove victim to fresh air and keep at rest in a position comfortable for breathing (P304+P340).

**If in eyes:** Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do so. Continue rinsing (P305+P351+P338).

If on skin (or hair): Immediately remove all contaminated clothing. Rinse skin with water (P303+P361+P353).

**If swallowed:** Rinse mouth. Do NOT induce vomiting (P301+P330+P331).

#### **SECTION 5 — FIRE FIGHTING MEASURES**

Nonflammable, noncombustible solution.

In case of fire: Use a tri-class dry chemical fire extinguisher.

NFPA CODE

H - 3 F - 0

D 0

R - 0

#### **SECTION 6 — ACCIDENTAL RELEASE MEASURES**

Remove all ignition sources and ventilate area. Restrict unprotected personnel from the area. Contain the spill with sand or absorbent material and deposit in a sealed bag or container. See Sections 8 and 13 for further information.

**Safety Data Sheet** 

Mercury(I) Nitrate Solution

Revision Date: March 25, 2014

**SDS #:** 505.10

#### **SECTION 7 — HANDLING AND STORAGE**

Flinn Suggested Chemical Storage Pattern: Inorganic #3. Store with amides, nitrates, nitrites and azides.

Store in a poisons cabinet. Never store with acetic acid. Keep and store away from clothing and combustible materials (P220). Take any precautions to avoid mixing with combustibles (P221). Use only in a hood or well-ventilated area (P271).

#### SECTION 8 — EXPOSURE CONTROLS, PERSONAL PROTECTION

Wear protective gloves, protective clothing, and eye protection (P280). Wash hands thoroughly after handling (P264). Use only in a hood or well-ventilated area (P271).

Exposure guidelines: (as mercury nitrate) PEL 0.1 mg/m<sup>3</sup> (OSHA); TLV 0.025 mg/m<sup>3</sup> (ACGIH) Readily absorbed through the skin.

#### **SECTION 9 — PHYSICAL AND CHEMICAL PROPERTIES**

Colorless liquid. Odorless.

#### **SECTION 10 — STABILITY AND REACTIVITY**

Avoid contact with organic materials. Shelf life: Indefinite, if stored properly.

#### **SECTION 11 — TOXICOLOGICAL INFORMATION**

Acute effects: Highly toxic, irritant, sensitizer ORL-RAT LD<sub>50</sub>: 170 mg/kg as mercury nitrate

Chronic effects: N.A. IHL-RAT  $LC_{50}$ : N.A. Target organs: N.A. SKN-RBT  $LD_{50}$ : N.A.

N.A. = Not available, not all health aspects of this substance have been fully investigated.

#### **SECTION 12 — ECOLOGICAL INFORMATION**

Very harmful to aquatic life. Avoid Release into the environment.

### <u>SECTION 13 — DISPOSAL CONSIDERATIONS</u>

Please review all federal, state and local regulations that may apply before proceeding.

Flinn Suggested Disposal Method #27b is one option.

#### **SECTION 14 — TRANSPORT INFORMATION**

Shipping name: Mercurous nitrate; Hazard class: 6.1 Poison; UN number: UN1627

N/A = Not applicable

#### **SECTION 15 — REGULATORY INFORMATION**

Not listed.

#### **SECTION 16 — OTHER INFORMATION**

This Safety Data Sheet (SDS) is for guidance and is based upon information and tests believed to be reliable. Flinn Scientific, Inc. makes no guarantee of the accuracy or completeness of the data and shall not be liable for any damages relating thereto. The data is offered solely for your consideration, investigation, and verification. The data should not be confused with local, state, federal or insurance mandates, regulations, or requirements and CONSTITUTE NO WARRANTY. Any use of this data and information must be determined by the science instructor to be in accordance with applicable local, state or federal laws and regulations. The conditions or methods of handling, storage, use and disposal of the product(s) described are beyond the control of Flinn Scientific, Inc. and may be beyond our knowledge. FOR THIS AND OTHER REASONS, WE DO NOT ASSUME RESPONSIBILITY AND EXPRESSLY DISCLAIM LIABILITY FOR LOSS, DAMAGE OR EXPENSE ARISING OUT OF OR IN ANY WAY CONNECTED WITH THE HANDLING, STORAGE, USE OR DISPOSAL OF THIS PRODUCT(S).

Consult your copy of the Flinn Science Catalog/Reference Manual for additional information about laboratory chemicals.

Revision Date: March 25, 2014